ClearWayTM RX

Rapid Exchange Therapeutic Perfusion Catheter

Instructions For Use

ATRIUM

www.atriummed.com

DESCRIPTION

The Atrium ClearWay™ RX is a therapeutic perfusion catheter for localized delivery of therapeutic agents to the coronary or peripheral vasculature. This rapid exchange catheter is .014" guidewire compatible. The 24cm distal portion of the catheter is dual lumen with the outer lumen for inflation of the balloon and inner lumen as a guidewire channel. ClearWay RX perfuses fluid out through the microporous PTFE material during continuous low-pressure injection (1-4 ATM). Radiopaque marker(s) are located within the balloon treatment zone to assist in catheter placement under fluoroscopy. This catheter is non-pyrogenic and is for single patient use only.

INDICATIONS

The ClearWay RX is indicated for localized perfusion of various diagnostic and therapeutic agents into the coronary and peripheral vasculature. It is not indicated for use in the neurovasculature.

CONTRAINDICATIONS

Do not use this device for percutaneous transluminal angioplasty or stent deployment.

WARNINGS

1. The microporous PTFE material must be purged of all air prior to use with the distal protective sheath in place.
2. To avoid air emboli, do not use air to inflate the microporous PTFE material.
3. To minimize the risk of vessel damage or balloon damage, do not exceed the maximum recommended inflation pressure provided on the package label.
4. All agents used with this device must be used in accordance with the manufacturer's instructions for use.
5. Do not expose this catheter to organic solvents or ionizing radiation.
6. Remove the distal protective sheath prior to catheter insertion.

PRECAUTIONS

1. Only physicians trained in surgical, angiographic and percutaneous interventional techniques should use this device.
2. The enclosed instructions for use provide technical guidance only, and do not eliminate the need for formal training in the use of interventional devices.
3. Do not force catheter passage or withdrawal without first determining the cause of catheter resistance under fluoroscopy.
4. Do not use this device if the sterile package has become damaged or previously opened.
5. Do not use automatic or high pressure injection systems as they may exceed maximum recommended inflation pressures and damage the PTFE material, catheter shaft or both.
6. ClearWay RX is intended for single site delivery. Do not reposition for reuse at another site.
7. Once withdrawn from the body, do not reintroduce the ClearWay RX catheter.

POSSIBLE COMPLICATIONS

- Hematoma at the access point
- Acute myocardial infarction
- Spontaneous thrombus formation and occlusion of the artery, vein or graft
- Cerebral circulatory disorders
- Unknown side effects with regard to various irrigating solutions and concentrations
- Possible flow surface disruption
- Hypotension
- Ischemia
- Palpitations
- Vessel dissection
- Anastomotic rupture or injury
- Restenosis
- Hemorrhage or hematoma
- Unstable angina
- Arrhythmias
- Allergic reaction to contrast medium

- Infection
- Coronary artery spasm
- Embolism
- Death

DIRECTIONS FOR USE

A. BALLOON PREPARATION

1. Remove catheter from the package.
2. A 10cc or larger syringe, or low pressure inflation device is recommended for use with this catheter.
3. The ClearWay RX catheter must be primed and purged of all air with sterile fluid prior to patient use with the distal protective sheath in place. Do not exceed 1 ATM.
4. The distal protective sheath must be kept in place over the PTFE balloon material during all purging steps.
5. Do not pull a vacuum on the syringe when purging the PTFE balloon material.
6. Prior to insertion remove the distal protective sheath and visually inspect the PTFE balloon material for any damage.


B. CATHETER INSERTION

1. See corresponding package label for guidewire and guide catheter compatibility.
2. If necessary, the guidewire lumen can be irrigated with sterile saline solution prior to patient use.
3. Carefully introduce the ClearWay RX catheter tip onto the guidewire and insert.
4. Under fluoroscopy advance the microporous balloon segment to the desired treatment location. Once placement is confirmed using the radiopaque markers, infusion can begin.
5. Infuse prescribed amount of fluid. Do not exceed the maximum recommended pressure found on the corresponding package label.
6. Once the prescribed volume of fluid has been infused, immediately release syringe pressure by pulling back on the syringe barrel.
7. After infusion is complete, ensure that the device is deflated prior to withdrawal.
8. The ClearWay RX catheter is for single site infusion only.

STORAGE

Store in a cool, dry place away from intense light.

SYMBOLS USED ON PRODUCT PACKAGING

©2007 All Rights Reserved.

This device is covered under the following U.S. patent: 6395208B1

Atrium and ClearWay are trademarks of Atrium Medical Corporation.

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©2007 All rights reserved.
Hudson, New Hampshire 03051 U.S.A.
ATRIUM EUROPE B.V.
ATRIUM AUSTRALIA-PACIFIC RIM PTY. LTD.